

PVI00590EN_V01

Endless Foil System

For closed filling of flexible containers

- Endless foil system for dust-free filling of pharmacological and toxic, hazardous materials
- High containment level to OEB 5 (OEL < 1 µg/m³)
- Highest operator and product protection
- Simple to operate, reliable handling and safe to use
- Endless foil system, closure system and filling equipment from one source
- Ex-version according to ATEX available
- Anti-static foil packages available
- Different qualities of foil packages available: e.g. clean room packed or not.



Function

The product is fed into the bag by means of a metering device bag until the necessary volume or weight is reached. The foil is then pulled down a little and then closed 2 times by means of welding scissors, pneumatic clippers, band rappers, Novi binders or similar. (or once with the cutting/welding scissors) and separated between the seals.

For the completely closed liner change, the liner carrier is equipped with the proven safe-change system.

